

Product datasheet for **HP205100**

APBA1 Human qPCR Primer Pair (NM_001163)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	320
Forward Sequence:	GGCGAATCTCAGCCAATGACTG
Reverse Sequence:	GTCCGCAATGTAGGAAATGGTCC
ACCN:	<u>NM_001163</u> , <u>NM_001163.1</u> , <u>NM_001163.2</u> , <u>NM_001163.3</u> , <u>BC156047</u> , <u>BC157030</u> , <u>BF435835</u> , <u>BQ269294</u> , <u>NM_001163.4</u>
UniProt ID:	<u>Q02410</u>
Synonyms:	D9S411E; LIN10; MINT1; X11; X11A; X11ALPHA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)