

## Product datasheet for **HP205101**

### FE65 (APBB1) Human qPCR Primer Pair (NM\_001164)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	322
Forward Sequence:	CGCTGCTTGGTAAATGGACTCTC
Reverse Sequence:	GGACCAACTCATTCTTAGCGC
ACCN:	<a href="#">NM_001164</a> , <a href="#">NM_001164.1</a> , <a href="#">NM_001164.2</a> , <a href="#">NM_001164.3</a> , <a href="#">NM_001164.4</a> , <a href="#">BC010854</a> , <a href="#">BG333769</a> , <a href="#">BG699532</a> , <a href="#">BM145723</a> , <a href="#">BX538185</a> , <a href="#">NM_001164.5</a>
UniProt ID:	<a href="#">O00213</a>
Synonyms:	FE65; MGC:9072; RIR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)