

Product datasheet for **HP205117**

Zinc Alpha 2 Glycoprotein (AZGP1) Human qPCR Primer Pair (NM_001185)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	563
Forward Sequence:	AGGACAGCCAACCTTCAGAAGGC
Reverse Sequence:	CGATCTCACAACCAAACCTTCCC
ACCN:	NM_001185 , NM_001185.1 , NM_001185.2 , NM_001185.3 , BC033830 , BC033830.1 , BC005306 , BC005908 , BC014470 , BM796305
UniProt ID:	P25311
Synonyms:	ZA2G; ZAG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)