

Product datasheet for **HP205121**

BCAT2 Human qPCR Primer Pair (NM_001190)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	587
Forward Sequence:	CGCTGAATGGTGTATCCTGCC
Reverse Sequence:	CAGCAACTGCTTCATGGTGATCG
ACCN:	NM_001190 , NM_001190.1 , NM_001190.2 , NM_001190.3 , BC001900 , BC001900.1 , BC004243 , BE139200 , BM802402 , BU178552 , NM_001190.4
UniProt ID:	O15382
Synonyms:	BCAM; BCATM; BCT2; HVLI; PP18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)