

Product datasheet for **HP205123**

BCMA (TNFRSF17) Human qPCR Primer Pair (NM_001192)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	608
Forward Sequence:	TCTTTGGCAGTTTTCGTGCTAATG
Reverse Sequence:	CCAGGTCAATGTTAGCCATGCC
ACCN:	<u>NM_001192</u> , <u>NM_001192.1</u> , <u>NM_001192.2</u> , <u>BC058291</u> , <u>BC058291.1</u> , <u>NM_001192.3</u>
UniProt ID:	<u>Q02223</u>
Synonyms:	BCM; BCMA; CD269; TNFRSF13A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)