

Product datasheet for **HP205127**

BIK Human qPCR Primer Pair (NM_001197)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	638
Forward Sequence:	GGAGGTTCTTGGCATGACTGAC
Reverse Sequence:	TGAGGCTCACGTCCATCTCGTC
ACCN:	<u>NM_001197</u> , <u>NM_001197.1</u> , <u>NM_001197.2</u> , <u>NM_001197.3</u> , <u>NM_001197.4</u> , <u>BC001599</u> , <u>BC001599.1</u> , <u>BM765910</u> , <u>BP329261</u> , <u>BT006728</u> , <u>NM_001197.5</u>
UniProt ID:	<u>Q13323</u>
Synonyms:	BIP1; BP4; NBK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)