

## Product datasheet for **HP205132**

### **BMP4 Human qPCR Primer Pair (NM\_001202)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	652
<b>Forward Sequence:</b>	CTGGTCTTGAGTATCCTGAGCG
<b>Reverse Sequence:</b>	TCACCTCGTTCTCAGGGATGCT
<b>ACCN:</b>	<a href="#">NM_001202</a> , <a href="#">NM_001202.1</a> , <a href="#">NM_001202.2</a> , <a href="#">NM_001202.3</a> , <a href="#">NM_001202.4</a> , <a href="#">NM_001202.5</a> , <a href="#">BC020546</a> , <a href="#">BC020546.2</a> , <a href="#">BE877424</a> , <a href="#">BQ083000</a> , <a href="#">BU683974</a> , <a href="#">BX161385</a> , <a href="#">BX161438</a>
<b>UniProt ID:</b>	<a href="#">P12644</a>
<b>Synonyms:</b>	BMP2B; BMP2B1; MCOPS6; OFC11; ZYME
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)