

Product datasheet for **HP205153**

Caspase-6 (CASP6) Human qPCR Primer Pair (NM_001226)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	839
Forward Sequence:	AGGTGGATGCAGCCTCCGTTTA
Reverse Sequence:	ATGAGCCGTTACAGTTTCCCG
ACCN:	NM_001226 , NM_001226.1 , NM_001226.2 , NM_001226.3 , BC004460 , BC004460.1 , BC000305 , BM837197 , BX380795 , NM_001226.4
UniProt ID:	P55212
Synonyms:	MCH2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)