

Product datasheet for **HP205159**

Caveolin 2 (CAV2) Human qPCR Primer Pair (NM_001233)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	858
Forward Sequence:	TTCTCTTTGCCACCCTCAGCTG
Reverse Sequence:	GAAGCATCGTCCTACGCTCGTA
ACCN:	NM_001233 , NM_001233.1 , NM_001233.2 , NM_001233.3 , NM_001233.4 , BC005256 , BC005256.1 , BM542127 , BM762545 , BT007051 , BX350116 , NM_001233.5
UniProt ID:	P51636
Synonyms:	CAV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)