

Product datasheet for **HP205169**

CD30 (TNFRSF8) Human qPCR Primer Pair (NM_001243)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	943
Forward Sequence:	ATCTGTGCCACATCAGCCACCA
Reverse Sequence:	AAGGTGGTGTCTCTCAGCCA
ACCN:	<u>NM_001243</u> , <u>NM_001243.1</u> , <u>NM_001243.2</u> , <u>NM_001243.3</u> , <u>NM_001243.4</u> , <u>BC063482</u> , <u>BC073839</u> , <u>BC136400</u> , <u>BC141804</u> , <u>BC171762</u> , <u>NM_001243.5</u>
UniProt ID:	<u>P28908</u>
Synonyms:	CD30; D1S166E; Ki-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)