

## Product datasheet for **HP205171**

### SIGLEC6 Human qPCR Primer Pair (NM\_001245)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	946
Forward Sequence:	TTTCACCTGCCGTGCTCAGCAT
Reverse Sequence:	ACCAGGGTTGTGATGCTAGCTC
ACCN:	<u><a href="#">NM_001245</a></u> , <u><a href="#">NM_001245.1</a></u> , <u><a href="#">NM_001245.2</a></u> , <u><a href="#">NM_001245.3</a></u> , <u><a href="#">NM_001245.4</a></u> , <u><a href="#">NM_001245.5</a></u> , <u><a href="#">NM_001245.6</a></u> , <u><a href="#">BC035359</a></u> , <u><a href="#">BC035359.1</a></u> , <u><a href="#">BC140798</a></u> , <u><a href="#">NM_001245.7</a></u>
UniProt ID:	<u><a href="#">O43699</a></u>
Synonyms:	CD33L; CD33L1; CD33L2; CD327; CDW327; OBBP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)