

## Product datasheet for **HP205180**

### CDHH (CDH13) Human qPCR Primer Pair (NM\_001257)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1012
Forward Sequence:	TGGCAGAACTCGTGATTGTCGG
Reverse Sequence:	GGGAGATACCACAATGGACCTC
ACCN:	<a href="#">NM_001257</a> , <a href="#">NM_001257.2</a> , <a href="#">NM_001257.3</a> , <a href="#">NM_001257.4</a> , <a href="#">BC030653</a> , <a href="#">BC012615</a> , <a href="#">BC028624</a> , <a href="#">BG911233</a> , <a href="#">BI756586</a> , <a href="#">BU162665</a> , <a href="#">BU627784</a> , <a href="#">NM_001257.5</a>
UniProt ID:	<a href="#">P55290</a>
Synonyms:	CDHH; P105
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)