

## Product datasheet for **HP205183**

### CDK8 Human qPCR Primer Pair (NM\_001260)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1024
Forward Sequence:	GCTGATAGGAAGGTGTGGCTTC
Reverse Sequence:	CCGAGGTAAGTGAAGTGGCTTC
ACCN:	<a href="#">NM_001260</a> , <a href="#">NM_001260.1</a> , <a href="#">NM_001260.2</a> , <a href="#">BC069634</a> , <a href="#">BC094784</a> , <a href="#">BC104492</a> , <a href="#">BC105950</a> , <a href="#">BC107601</a> , <a href="#">BM021447</a> , <a href="#">BM021938</a> , <a href="#">NM_001260.3</a>
UniProt ID:	<a href="#">P49336</a>
Synonyms:	IDDHBA; K35
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)