

Product datasheet for **HP205184**

CDK9 Human qPCR Primer Pair (NM_001261)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1025
Forward Sequence:	CCATTACAGCCTTGCGGGAGAT
Reverse Sequence:	CAGCAAGGTCATGCTCGCAGAA
ACCN:	<u>NM_001261</u> , <u>NM_001261.1</u> , <u>NM_001261.2</u> , <u>NM_001261.3</u> , <u>BC001968</u> , <u>BC001968.1</u> , <u>BE549667</u> , <u>BT019903</u> , <u>BX091973</u> , <u>NM_001261.4</u>
UniProt ID:	<u>P50750</u>
Synonyms:	C-2k; CDC2L4; CTK1; PITALRE; TAK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)