

Product datasheet for **HP205193**

Chromogranin A (CHGA) Human qPCR Primer Pair (NM_001275)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1113
Forward Sequence:	GGTTCTTGAGAACCAGAGCAGC
Reverse Sequence:	GCTTCACCACTTTTCTCTGCCTC
ACCN:	NM_001275 , NM_001275.1 , NM_001275.2 , NM_001275.3 , BC006459 , BC006459.1 , BC001059 , BC009384 , BC012755 , BT006869 , BX248263 , NM_001275.4
UniProt ID:	P10645
Synonyms:	CGA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)