

Product datasheet for **HP205198**

CIRBP Human qPCR Primer Pair (NM_001280)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1153
Forward Sequence:	GAGTCAGAGTGGTGGCTACAGT
Reverse Sequence:	CACAGCCATTGGAAGGACGATC
ACCN:	<u>NM_001280</u> , <u>NM_001280.1</u> , <u>NM_001280.2</u> , <u>BC000901</u> , <u>BC000901.1</u> , <u>BC000403</u> , <u>BI819825</u> , <u>BM706055</u> , <u>BP195090</u> , <u>BT007210</u> , <u>BU689520</u>
UniProt ID:	<u>Q14011</u>
Synonyms:	CIRP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)