

Product datasheet for **HP205209**

CLPTM1 Human qPCR Primer Pair (NM_001294)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1209
Forward Sequence:	GCTCACCATCATCGTGTCTATCG
Reverse Sequence:	CCACGAATGACTGGAAAACGCC
ACCN:	NM_001294 , NM_001294.1 , NM_001294.2 , NM_001294.3 , BC012359 , BC012359.1 , BC004865 , BI916245 , BM676309 , BQ574654 , BT007262 , NM_001294.4
UniProt ID:	O96005
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)