

Product datasheet for **HP205210**

CCR1 Human qPCR Primer Pair (NM_001295)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1230
Forward Sequence:	CAACTCCGTGCCAGAAGGTGAA
Reverse Sequence:	G TTCAGGAGGTAGATGCTGGTC
ACCN:	<u>NM_001295</u> , <u>NM_001295.1</u> , <u>NM_001295.2</u> , <u>BC051306</u> , <u>BC051306.1</u> , <u>BC064991</u> , <u>BT019908</u> , <u>BT019909</u> , <u>BX337412</u> , <u>NM_001295.3</u>
UniProt ID:	<u>P32246</u>
Synonyms:	CD191; CKR-1; CKR1; CMKBR1; HM145; MIP1aR; SCYAR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)