

Product datasheet for **HP205213**

Calponin (CNN1) Human qPCR Primer Pair (NM_001299)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1264
Forward Sequence:	CCAACGACCTGTTTGAGAACACC
Reverse Sequence:	ATTTCGCTCCTGCTTCTCTGC
ACCN:	<u>NM_001299</u> , <u>NM_001299.1</u> , <u>NM_001299.2</u> , <u>NM_001299.3</u> , <u>NM_001299.4</u> , <u>NM_001299.5</u> , <u>BC022015</u> , <u>BC022015.2</u> , <u>BC036307</u> , <u>BG718473</u> , <u>NM_001299.6</u>
UniProt ID:	<u>P51911</u>
Synonyms:	HEL-S-14; Sm-Calp; SMCC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)