

## Product datasheet for **HP205217**

### Claudin 4 (CLDN4) Human qPCR Primer Pair (NM\_001305)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1364
Forward Sequence:	AGTGCAAGGTGTACGACTCGCT
Reverse Sequence:	CGCTTTCATCCTCCAGGCAGTT
ACCN:	<a href="#">NM_001305</a> , <a href="#">NM_001305.1</a> , <a href="#">NM_001305.2</a> , <a href="#">NM_001305.3</a> , <a href="#">NM_001305.4</a> , <a href="#">BC000671</a> , <a href="#">BC000671.2</a> , <a href="#">BM768799</a> , <a href="#">BT006989</a> , <a href="#">NM_001305.5</a>
UniProt ID:	<a href="#">O14493</a>
Synonyms:	CPE-R; CPER; CPETR; CPETR1; hCPE-R; WBSCR8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)