

Product datasheet for **HP205222**

CRIP1 Human qPCR Primer Pair (NM_001311)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1396
Forward Sequence:	GCCTGAAGTGCAGAAATGTGG
Reverse Sequence:	CCCAAACATGGCTGCGTAGCAG
ACCN:	<u>NM_001311</u> , <u>NM_001311.1</u> , <u>NM_001311.2</u> , <u>NM_001311.3</u> , <u>NM_001311.4</u> , <u>BC002738</u> , <u>BC002738.2</u> , <u>BQ128075</u> , <u>BT007181</u> , <u>BU500642</u> , <u>NM_001311.5</u>
UniProt ID:	<u>P50238</u>
Synonyms:	CRHP; CRIP; CRP-1; CRP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)