

Product datasheet for **HP205228**

CSNK1G2 Human qPCR Primer Pair (NM_001319)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1455
Forward Sequence:	AAGGAGCGGTACCAGAAGATCG
Reverse Sequence:	GAAGAGCTCCGCAGGTAGTCA
ACCN:	<u>NM_001319</u> , <u>NM_001319.1</u> , <u>NM_001319.2</u> , <u>NM_001319.4</u> , <u>NM_001319.5</u> , <u>NM_001319.6</u> , <u>BC020972</u> , <u>BC020972.1</u> , <u>BC018693</u> , <u>BC018699</u> , <u>BT009922</u>
UniProt ID:	<u>P78368</u>
Synonyms:	CK1g2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)