

Product datasheet for **HP205239**

Cathepsin W (CTSW) Human qPCR Primer Pair (NM_001335)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1521
Forward Sequence:	CGTGACCATCAACATGAAGCCC
Reverse Sequence:	CCTCTGACTTGACGCTGCCAAA
ACCN:	<u>BC048255</u> , <u>NM_001335</u> , <u>NM_001335.2</u> , <u>NM_001335.3</u> , <u>BC048255.1</u> , <u>BC035637</u> , <u>BI833187</u> , <u>BM996548</u>
UniProt ID:	<u>P56202</u>
Synonyms:	LYPN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)