

Product datasheet for **HP205275**

OVCA1 (DPH1) Human qPCR Primer Pair (NM_001383)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1801
Forward Sequence:	CACTACACACCTCCTGGACTCT
Reverse Sequence:	CTCACACGATACTCGGCTTTCAG
ACCN:	NM_001383 , NM_001383.1 , NM_001383.2 , NM_001383.3 , NM_001383.4 , BC003099 , BC013205 , BC096086 , BC096087 , BC096088 , BC096089 , BP417755 , BT019878
UniProt ID:	Q9BZG8
Synonyms:	DEDSSH; DPH2L; DPH2L1; OVCA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)