

Product datasheet for **HP205281**

Dystrobrevin alpha (DTNA) Human qPCR Primer Pair (NM_001390)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1837
Forward Sequence:	CGGCTTGATGAAGAACACAGGC
Reverse Sequence:	GATGTCAGGAGCACTTCTCTGC
ACCN:	NM_001390 , NM_001390.1 , NM_001390.2 , NM_001390.3 , NM_001390.4 , BC005300 , BC033969 , BC047095 , BG250424 , BG705393 , BI668892 , BT006937 , BT007812 , BX103860
UniProt ID:	Q9Y4J8
Synonyms:	D18S892E; DRP3; DTN; DTN-A; LVNC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)