

Product datasheet for **HP205287**

EDG1 (S1PR1) Human qPCR Primer Pair (NM_001400)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1901
Forward Sequence:	CCTGTGACATCCTCTTCAGAGC
Reverse Sequence:	CACTTGCAGCAGGACATGATCC
ACCN:	<u>NM_001400</u> , <u>NM_001400.1</u> , <u>NM_001400.2</u> , <u>NM_001400.3</u> , <u>NM_001400.4</u> , <u>BC018650</u> , <u>BC018650.1</u> , <u>BI598518</u> , <u>BP236900</u> , <u>BU677564</u> , <u>NM_001400.5</u>
UniProt ID:	<u>P21453</u>
Synonyms:	CD363; CHEDG1; D1S3362; ECGF1; EDG-1; EDG1; S1P1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)