

Product datasheet for **HP205288**

EDG2 (LPAR1) Human qPCR Primer Pair (NM_001401)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1902
Forward Sequence:	GGCTATGTTGCCAGAGGACTA
Reverse Sequence:	GGAGTCCAGCAGATGATAAAGGC
ACCN:	<u>NM_001401</u> , <u>NM_001401.1</u> , <u>NM_001401.2</u> , <u>NM_001401.3</u> , <u>BC030615</u> , <u>BC030615.2</u> , <u>BC036034</u> , <u>BC036034.1</u> , <u>BC073167</u> , <u>BI461641</u> , <u>NM_001401.5</u>
UniProt ID:	<u>Q92633</u>
Synonyms:	edg-2; EDG2; Gpcr26; LPA1; Mrec1.3; rec.1.3; vzg-1; VZG1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)