

Product datasheet for **HP205302**

ELAVL1 Human qPCR Primer Pair (NM_001419)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1994
Forward Sequence:	TGTTCTCTCGGTTTGGCGGAT
Reverse Sequence:	TCTTCTGCCTCCGACCGTTGT
ACCN:	NM_001419 , NM_001419.1 , NM_001419.2 , BC003376 , BC003376.2 , BM664469 , BM696191 , BT009793 , BU542524 , BX647856 , NM_001419.3
UniProt ID:	Q15717
Synonyms:	ELAV1; Hua; HUR; MeIG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)