

Product datasheet for **HP205313**

HIF 2 alpha (EPAS1) Human qPCR Primer Pair (NM_001430)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2034
Forward Sequence:	CTGTGTCTGAGAAGAGTAACTTCC
Reverse Sequence:	TTGCCATAGGCTGAGGACTCCT
ACCN:	<u>NM_001430</u> , <u>NM_001430.1</u> , <u>NM_001430.2</u> , <u>NM_001430.3</u> , <u>NM_001430.4</u> , <u>BC051338</u> , <u>BC051338.1</u> , <u>BC015869</u> , <u>BM671511</u> , <u>NM_001430.5</u>
UniProt ID:	<u>Q99814</u>
Synonyms:	bHLHe73; ECT4; HIF2A; HLF; MOP2; PASD2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)