

## Product datasheet for **HP205314**

### EPB4IL2 (EPB41L2) Human qPCR Primer Pair (NM\_001431)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                              |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                            |
| Gene ID:          | 2037                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | GTAAACGGGTCTCCAGGAGTCT                                                                                                                                                                                                                                                                                                                       |
| Reverse Sequence: | ACCACGGCAATGCTGACAAGTC                                                                                                                                                                                                                                                                                                                       |
| ACCN:             | <u><a href="#">NM_001431</a></u> , <u><a href="#">NM_001431.1</a></u> , <u><a href="#">NM_001431.2</a></u> , <u><a href="#">NM_001431.3</a></u> , <u><a href="#">BC034718</a></u> , <u><a href="#">BI461477</a></u> , <u><a href="#">BM679819</a></u> , <u><a href="#">BQ446747</a></u> , <u><a href="#">NM_001431.4</a></u>                 |
| UniProt ID:       | <u><a href="#">O43491</a></u>                                                                                                                                                                                                                                                                                                                |
| Synonyms:         | 4.1-G; 4.1G                                                                                                                                                                                                                                                                                                                                  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)                                                                                                                                                                                                                                               |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                              |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                 |



[View online »](#)