

Product datasheet for **HP205317**

Fibrillarin (FBL) Human qPCR Primer Pair (NM_001436)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2091
Forward Sequence:	GAGGCTTCCATTCTGGTGGCAA
Reverse Sequence:	CAGGTTCTGGTGACCAAGTGCA
ACCN:	NM_001436 , NM_001436.1 , NM_001436.2 , NM_001436.3 , BC019260 , BC019260.1 , BC015218 , BC019609 , BG831506 , BT006830 , BT020144 , NM_001436.4
UniProt ID:	P22087
Synonyms:	FIB; FLRN; Nop1; RNU3IP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)