

## Product datasheet for **HP205318**

### Estrogen Related Receptor gamma (ESRRG) Human qPCR Primer Pair (NM\_001438)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2104
Forward Sequence:	CGCAGGATAGATGCGGAGAACA
Reverse Sequence:	TTCAGCCACCAACAAATGTGAGAC
ACCN:	<a href="#">NM_001438</a> , <a href="#">NM_001438.1</a> , <a href="#">NM_001438.2</a> , <a href="#">NM_001438.3</a> , <a href="#">BC008218</a> , <a href="#">BC064700</a> , <a href="#">BM723962</a> , <a href="#">BP349857</a> , <a href="#">BU752915</a> , <a href="#">BX486063</a> , <a href="#">BX647076</a> , <a href="#">BX649009</a> , <a href="#">NM_001438.4</a>
UniProt ID:	<a href="#">P62508</a>
Synonyms:	ERR-gamma; ERR3; ERRg; ERRgamma; NR3B3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)