

Product datasheet for **HP205332**

FOXO3 Human qPCR Primer Pair (NM_001455)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2309
Forward Sequence:	TCTACGAGTGGATGGTGC GTT G
Reverse Sequence:	CTCTTGCCAGTTCCTCATTCTG
ACCN:	<u>NM_001455</u> , <u>NM_001455.1</u> , <u>NM_001455.2</u> , <u>NM_001455.3</u> , <u>BC021224</u> , <u>BC021224.2</u> , <u>BC068552</u> , <u>BC020227</u> , <u>BC045800</u>
UniProt ID:	<u>O43524</u>
Synonyms:	AF6q21; FKHRL1; FKHRL1P2; FOXO2; FOXO3A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)