

Product datasheet for **HP205333**

Filamin B (FLNB) Human qPCR Primer Pair (NM_001457)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2317
Forward Sequence:	CCTTCAAGGTGGCTGCTCACTGA
Reverse Sequence:	CCCTCAACAGTTATGCCAAGCC
ACCN:	<u>NM_001457</u> , <u>NM_001457.1</u> , <u>NM_001457.2</u> , <u>NM_001457.3</u> , <u>BC156184</u> , <u>BX641085</u> , <u>NM_001457.4</u>
UniProt ID:	<u>O75369</u>
Synonyms:	ABP-278; ABP-280; AOI; FH1; FLN-B; FLN1L; LRS1; SCT; TABP; TAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)