

Product datasheet for **HP205334**

Filamin C (FLNC) Human qPCR Primer Pair (NM_001458)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2318
Forward Sequence:	ATGGTAGCTGCACCGTGGAGTA
Reverse Sequence:	TCCACCACATCCTTCACTGGCA
ACCN:	<u>NM_001458</u> , <u>NM_001458.1</u> , <u>NM_001458.2</u> , <u>NM_001458.3</u> , <u>NM_001458.4</u>
UniProt ID:	<u>Q14315</u>
Synonyms:	ABP-280; ABP280A; ABPA; ABPL; CMH26; FLN2; MFM5; MPD4; RCM5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)