

Product datasheet for **HP205335**

Flt3 ligand (FLT3LG) Human qPCR Primer Pair (NM_001459)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2323
Forward Sequence:	AGACTGTCGCTGGGTCCAAGAT
Reverse Sequence:	GGAGATGTTGGTCTGGACGAAG
ACCN:	<u>NM_001459</u> , <u>NM_001459.1</u> , <u>NM_001459.2</u> , <u>NM_001459.3</u> , <u>BC126293</u> , <u>BC006331</u> , <u>BC011914</u> , <u>BC028001</u> , <u>BC136464</u> , <u>BC144129</u> , <u>NM_001459.4</u>
UniProt ID:	<u>P49771</u>
Synonyms:	FL; FLG3L; FLT3L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)