

Product datasheet for **HP205344**

GABA B Receptor 1 (GABBR1) Human qPCR Primer Pair (NM_001470)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2550
Forward Sequence:	CCTGAACAAGACATCTGGAGGAG
Reverse Sequence:	GCTGGCATCAAACACCACATGG
ACCN:	<u>NM_001470</u> , <u>NM_001470.1</u> , <u>NM_001470.2</u> , <u>NM_001470.3</u> , <u>BC050532</u> , <u>BC050532.1</u> , <u>BC023336</u> , <u>BC041332</u> , <u>BC042598</u> , <u>NM_001470.4</u>
UniProt ID:	<u>Q9UBS5</u>
Synonyms:	GABABR1; GABBR1-3; GB1; GPRC3A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)