

Product datasheet for **HP205361**

GDF1 Human qPCR Primer Pair (NM_001492)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2657
Forward Sequence:	GTCACCCTGCAACCGTGCCAC
Reverse Sequence:	AGGTCGAAGACGACTGTCCACT
ACCN:	<u>NM_001492</u> , <u>NM_001492.1</u> , <u>NM_001492.2</u> , <u>NM_001492.3</u> , <u>NM_001492.4</u> , <u>NM_001492.5</u> , <u>BC022450</u> , <u>BC022450.1</u>
UniProt ID:	<u>P27539</u>
Synonyms:	CERS1; CHTD6; DORV; DTGA3; LAG1; LASS1; RAI; UOG1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)