

Product datasheet for **HP205380**

GTF2H2 Human qPCR Primer Pair (NM_001515)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2966
Forward Sequence:	CATGCCTGGACATACAAGTCGAG
Reverse Sequence:	CAGTGCAAACGCGAACTTCTGC
ACCN:	<u>NM_001515</u> , <u>NM_001515.2</u> , <u>NM_001515.3</u> , <u>BC005345</u> , <u>BC140303</u> , <u>BC141603</u> , <u>BG283896</u> , <u>BM083743</u> , <u>BT006773</u>
UniProt ID:	<u>Q13888</u>
Synonyms:	BTF2; BTF2 p44; BTF2P44; p44; T-BTF2P44; TFIIH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)