

Product datasheet for **HP205393**

HIF-1 alpha (HIF1A) Human qPCR Primer Pair (NM_001530)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 3091 |
| Forward Sequence: | TATGAGCCAGAAGAAGCTTTTAGGC |
| Reverse Sequence: | CACCTCTTTTGGCAAGCATCCTG |
| ACCN: | <u>NM_001530</u> , <u>NM_001530.1</u> , <u>NM_001530.2</u> , <u>NM_001530.3</u> , <u>BC012527</u> , <u>BC012527.2</u> , <u>BT009776</u> , <u>BX648795</u> , <u>NM_001530.4</u> |
| UniProt ID: | <u>Q16665</u> |
| Synonyms: | bHLHe78; HIF-1-alpha; HIF-1A; HIF-1alpha; HIF1; HIF1-ALPHA; MOP1; PASD8 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)