

Product datasheet for **HP205400**

HDJ2 (DNAJA1) Human qPCR Primer Pair (NM_001539)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3301
Forward Sequence:	GGTGAAGGAGACCAAGAACCAG
Reverse Sequence:	AGCCACACAGTGCTTCAACGAG
ACCN:	<u>NM_001539</u> , <u>NM_001539.1</u> , <u>NM_001539.2</u> , <u>NM_001539.3</u> , <u>BC008182</u> , <u>BC008182.1</u> , <u>BI462555</u> , <u>BM469197</u> , <u>BT007292</u> , <u>BU508753</u>
UniProt ID:	<u>P31689</u>
Synonyms:	DJ-2; DjA1; hDJ-2; HDJ2; HSDJ; HSJ-2; HSJ2; HSPF4; NEDD7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)