

## Product datasheet for **HP205406**

### IFIT2 Human qPCR Primer Pair (NM\_001547)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3433
Forward Sequence:	GGAGCAGATTCTGAGGCTTTGC
Reverse Sequence:	GGATGAGGCTTCCAGACTCCAA
ACCN:	<u><a href="#">NM_001547</a></u> , <u><a href="#">NM_001547.2</a></u> , <u><a href="#">NM_001547.3</a></u> , <u><a href="#">NM_001547.4</a></u> , <u><a href="#">BC032839</a></u> , <u><a href="#">BC032839.2</a></u> , <u><a href="#">BC005987</a></u>
UniProt ID:	<u><a href="#">P09913</a></u>
Synonyms:	cig42; G10P2; GARG-39; IFI-54; IFI-54K; IFI54; IFIT-2; ISG-54 K; ISG-54K; ISG54; P54
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)