

## Product datasheet for **HP205412**

### IKK beta (IKKB) Human qPCR Primer Pair (NM\_001556)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3551
Forward Sequence:	ACAGCGAGCAAACCGAGTTTGG
Reverse Sequence:	CCTCTGTAAGTCCACAATGTCGG
ACCN:	<a href="#">NM_001556</a> , <a href="#">NM_001556.1</a> , <a href="#">NM_001556.2</a> , <a href="#">BC108694</a> , <a href="#">BC006231</a> , <a href="#">BC090865</a> , <a href="#">BC127663</a> , <a href="#">BC141978</a> , <a href="#">BC144080</a> , <a href="#">BX648165</a> , <a href="#">NM_001556.3</a>
UniProt ID:	<a href="#">O14920</a>
Synonyms:	IKK-beta; IKK2; IKKB; IMD15; IMD15A; IMD15B; NFKBIKB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)