

Product datasheet for **HP205424**

EIF3E Human qPCR Primer Pair (NM_001568)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3646
Forward Sequence:	CTTCTGTGCGATCCACCACTG
Reverse Sequence:	TCAATCTGGCATCCAGTCTTGC
ACCN:	<u>NM_001568</u> , <u>NM_001568.1</u> , <u>NM_001568.2</u> , <u>BC017887</u> , <u>BC017887.1</u> , <u>BC000734</u> , <u>BC005944</u> , <u>BC008419</u> , <u>BC016706</u> , <u>BC021679</u> , <u>BG498300</u> , <u>NM_001568.3</u>
UniProt ID:	<u>P60228</u>
Synonyms:	eIF3-p46; EIF3-P48; EIF3S6; INT6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)