

Product datasheet for **HP205432**

ACAA1 Human qPCR Primer Pair (NM_001607)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	30
Forward Sequence:	GACAGGTCATCACGCTGCTCAA
Reverse Sequence:	CCAGGGTATTCAAAGACGGCAG
ACCN:	NM_001607 , NM_001607.1 , NM_001607.2 , NM_001607.3 , BC000635 , BC000635.2 , BC011977 , BC014474 , BC025780 , BC039837 , NM_001607.4
UniProt ID:	P09110
Synonyms:	ACAA; PTHIO; THIO
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)