

Product datasheet for **HP205434**

ACADSB Human qPCR Primer Pair (NM_001609)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	36
Forward Sequence:	GCCACCTATTTGCCTCAGCTCA
Reverse Sequence:	GTCAGCACTGCTGATCCACAT
ACCN:	NM_001609 , NM_001609.1 , NM_001609.2 , NM_001609.3 , BC013756 , BC013756.1 , BC020931 , BC032395 , BC045723 , BE875638 , BQ312626 , NM_001609.4
UniProt ID:	P45954
Synonyms:	2-MEBCAD; ACAD7; SBCAD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)