

Product datasheet for **HP205442**

PARP1 Human qPCR Primer Pair (NM_001618)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	142
Forward Sequence:	CCAAGCCAGTTCAGGACCTCAT
Reverse Sequence:	GGATCTGCCTTTTGCTCAGCTTC
ACCN:	<u>NM_001618</u> , <u>NM_001618.1</u> , <u>NM_001618.2</u> , <u>NM_001618.3</u> , <u>BC037545</u> , <u>BC037545.1</u> , <u>BC008660</u> , <u>BC014206</u> , <u>BC018620</u> , <u>BC021045</u>
UniProt ID:	<u>P09874</u>
Synonyms:	ADPRT; ADPRT 1; ADPRT1; ARTD1; pADPRT-1; PARP; PARP-1; PPOL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)