

Product datasheet for **HP205453**

FLAP (ALOX5AP) Human qPCR Primer Pair (NM_001629)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	241
Forward Sequence:	AAGTGGAGCACGAAAGCAGGAC
Reverse Sequence:	AGACCAGAGCACAGCGAGGAAA
ACCN:	NM_001629 , NM_001629.1 , NM_001629.2 , NM_001629.3 , BC018538 , BC018538.1 , BG569867 , BQ003388 , NM_001629.4
UniProt ID:	P20292
Synonyms:	FLAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)