

Product datasheet for **HP205464**

APE1 (APEX1) Human qPCR Primer Pair (NM_001641)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	328
Forward Sequence:	CTGCTCTTGAATGTGGATGGG
Reverse Sequence:	TCCAGGCAGCTCCTGAAGTTCA
Accession No:	<u>NM_001244249</u> , <u>NM_001641</u> , <u>NM_080648</u> , <u>NM_080649</u> , <u>BC004979</u> , <u>BC002338</u> , <u>BC095428</u> , <u>BC008145</u> , <u>BC019291</u> , <u>BT007236</u> , <u>BT020133</u> , <u>BX161454</u> , <u>NM_080648.3</u> , <u>NM_001244249.2</u> , <u>NM_080649.3</u> , <u>NM_001641.4</u>
Synonyms:	APE; APE1; APEN; APEX; APX; HAP1; REF1
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)