

## Product datasheet for **HP205468**

### Apolipoprotein CI (APOC1) Human qPCR Primer Pair (NM\_001645)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	341
Forward Sequence:	AGGACAAGGCTCGGGAACAT
Reverse Sequence:	GATGTCACCCTTCAGTCCTCA
ACCN:	<u><a href="#">NM_001645</a></u> , <u><a href="#">NM_001645.1</a></u> , <u><a href="#">NM_001645.2</a></u> , <u><a href="#">NM_001645.3</a></u> , <u><a href="#">NM_001645.4</a></u> , <u><a href="#">BC055093</a></u> , <u><a href="#">BC055093.1</a></u> , <u><a href="#">BC009698</a></u> , <u><a href="#">BT007142</a></u> , <u><a href="#">NM_001645.5</a></u>
UniProt ID:	<u><a href="#">P02654</a></u>
Synonyms:	Apo-CI; apo-CIB; ApoC-I; apoC-IB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)